



CALIBRATION INFORMATION

Things to Know:

- ➡ The Reagents Will Need to Be Calibrated if Using an Automated Testing System.
- ➡ Calibration solution is NOT SUPPLIED since the value of accuracy is dependent on what the user determines is the best standard of accuracy.

Calibration Options:

- 1 Users can perform manual tests of given aquarium water using their trusted test kits. Then calibrate using those values. Example home test kits include but are not limited to Salifert, Fauna Marin, Red Sea, API, etc.
- 2 Use calibration solutions from another preferred source/vendor. We recommend Fauna Marin's Multi-Reference Solution which is a ICP-tested reference solution with given values listed on the bottle.



With Fauna Marin Multi-Reference Solution, you can evaluate your home test kits under accurate conditions. The Reference Solution is based on a saltwater matrix and is quality tested using Fauna Marin's ICP-AES testing machine.

<https://cvue.us/multireference>

- 3 Use any water source with values you determine are accurate to your liking!



Need help...Visit :
support.coralvue.com